

## MEAS 400 and 700 Series Medical Probes

Cost-effective MEAS probes measure temperature with MEAS interchangeable thermistors, a highly accurate temperature technology that MEAS perfected over the decades. This thermal expertise, in combination with regulatory experience, manufacturing automation, and custom capability, positions MEAS as a leader in temperature measurement for healthcare facilities.

### MEAS 400 Series Reusable Temperature Probes


As demands on healthcare providers increase, healthcare professionals are taking the lead in finding answers to cost versus quality issues. MEAS 400 Series reusable probes are a sensible solution, offering substantial reductions in cost with no trade-off in quality.

### MEAS 400AC Series Medical Autoclavable Probes

MEAS 400AC Series Reusable Temperature Probes have the ability to withstand autoclave sterilization, generating less waste. Featuring gold-plated plugs and sealed thermistors, MEAS 400AC probes are resistant to standard chemicals and disinfectants. Guaranteed for 100 autoclave cycles or six months.

### MEAS 700 Series Reusable Probes

MEAS makes a complete line of reusable probes for 700 Series inputs providing linear response. Most MEAS 400 Series probes have a corresponding 700 Series probe. These are easily distinguished from the 400 Series probes at the connector end.



MEAS4012 Mini-Tempheart Shields



#### Dependability

- MEAS temperature probes deliver continuous, high-precision measurement

#### Repeatability

- Every MEAS probe is traceable to the U.S. National Institute of Standards and Technology (NIST)

#### Interchangeability

- MEAS patented manufacturing techniques ensure that probes are interchangeable within the same series



## MEAS 400 and 700 Series Probes Characteristics

<b>Range of Operation</b>	0° to 60°C
<b>Accuracy</b>	400 Series: ±0.2°C from -1° to 60°C, ±0.1°C from 25° to 45°C 700 Series: ±0.2°C from -1° to 60°C, ±0.15°C from 25° to 45°C
<b>Cleaning</b>	Probes should be cleaned with a mild detergent and water to remove excess bioburden and improve effectiveness of disinfection and sterilization.
<b>Disinfection</b>	Low level: cidex/glutaraldehyde High level: cidex/glutaraldehyde, dilute bleach, 70% isopropyl alcohol
<b>Sterilization</b>	Ethylene oxide gas, STERIS® System 1 MEAS 400AC only: steam autoclave (20 min. at 121-134 C), recommended ANSI/AAMI ST35
<b>Termination</b>	Standard .25" phone plug, compatible with most patient monitors
<b>Regulatory</b>	Compliant with FDA regulations, 93/42/EEC CE-Medical Device Directives, and EN12470

## MEAS 400 and 700 Series Probes Specifications

<b>400 Series</b>	401	Esophageal or Rectal (adult)	12 Fr.*, 4.0 mm tip diameter
	402	Esophageal or Rectal (pediatric)	10 Fr., 3.3 mm tip diameter
	403	Rectal (adult)	12 Fr., 4.0 mm diameter
	405	Gas/Air Probe	Measures or monitors ambient air temperature
	406	Rectal (pediatric)	9 Fr., 3.0 mm diameter
	409B	Skin (adult, flexible)	9.5 mm sensor disk diameter
	423	Rectal (pediatric, flexible)	10 Fr., probe length is 6.3 cm, diameter is 3.2 mm
	427	Skin (pediatric)	4.8 mm sensor disk diameter
<b>Autoclavable</b>	401AC	Esophageal or Rectal (adult)	14 Fr., 4.7 mm diameter
	402AC	Esophageal or Rectal (pediatric)	10 Fr., 3.3 mm diameter
	409AC	Skin (adult)	9.5 mm diameter disk
<b>700 Series</b>	701	Esophageal or Rectal (adult)	17 Fr., 6.3 mm tip diameter
	702A	Esophageal or Rectal (adult)	12 Fr., 4.0 mm tip diameter
	703	Rectal (adult)	12 Fr., 4.0 mm diameter
	708	Skin (adult)	9.5 mm sensor disk diameter
	709B	Skin (adult)	9.5 mm sensor disk diameter
	729	Skin (pediatric)	4.8 mm sensor disk diameter

\*Fr. = French Catheter scale

## Probe Accessories

<b>4009</b>	Tempheart Shields
<b>4012</b>	Mini-Tempheart Shields
Heat reflecting adhesive pads hold probes securely to the skin for accurate temperature readings	
<b>4930</b>	Adapter for Phillips Monitors

## Compatibility

- MEAS probes are designed to plug directly into standard MEAS 400 and 700 Series monitoring instruments

## Safety

- With no need for adapter boxes and electrical isolation of conductors, non-latex MEAS reusable probes virtually eliminate concerns about safety

## Cost-Effective

- MEAS reusable probes have a cost-per-use far below that of any competitive product or technology

To order or for more information, contact your local representative or the MEAS Temperature customer service team.

**937 427-1231**

800 747-5367 U.S. only  
937 427 1640 fax

[www.meas-spec.com](http://www.meas-spec.com)

2670 Indian Ripple Road  
Dayton, Ohio  
45440-3605 USA

BetaTHERM Sensors  
Ballybrit Business Park  
Galway, Ireland

ISO 9001  
ISO 13485



MEAS 400 Series Medical Probes